



#### NYLON PLASTIC COMPONENTS

##### For cleaning

Use soft cloth and non-abrasive neutral detergent. Do not use alcohol, solvent-based, acid or abrasive cleaners.

##### For disinfection

We recommend the use of sodium hypochlorite (0.1% - 0.5% sol.) or hydrogen peroxide.\*

**Before disinfecting the entire surface of any component, it is advisable to proceed with tests on a small part of the product.**

#### PRODUCT NAME

**Epsilon Universal plate**

#### TYPE

**Weight activated mechanism**

#### MATERIALS

**Nylon plastic**

\* Analysis of the effects resulting from the application of chemicals on product surfaces, performed by an accredited laboratory.